## TG - CUPTOARE TUBULARE CU GRADIENT

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The TG2 tube furnace includes a 25 mm long unheated zone barrier between the two 200 mm heated zones while the TG3 includes two 75 mm long unheated zone barriers between the three 150 mm heated zones. Each heated zone has its own temperature controller and thermocouple.

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HIGH-QUALITY THERMAL INSULATION


## LOW ENERGY CONSUMPTION

Carbolite's new tube furnaces are very energy-efficient. Their high-quality thermal insulation is designed for longevity, low case temperature and most of all outstandingly low energy consumption.

Using a Carbolite tube furnace reduces your energy cost as well as the carbon footprint.

**UP TO 50% LESS ENERGY REQUIRED!**

**COMPARISON WITH EQUIVALENT COMPETITOR PRODUCT**

Energy consumption for exemplary Carbolite TS1 12/60/600 model (red) and equivalent competitor product (yellow); measured at each temperature after 2 hours soak.

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## TG GRADIENT TUBE FURNACE RANGE - FAQ

### **WHAT IS A GRADIENT TUBE FURNACE?**

A gradient tube furnace is a split tube furnace with multiple heated zones, fitted with zone barriers. The zone barriers are pieces of insulation that fit between each heating element module, and help reduce the amount of heat transfer between heated zones.

### **HOW ARE THE GRADIENTS IN A GRADIENT TUBE FURNACE CONTROLLED?**

Each heated zone is controlled by its own independent temperature controller, so it is possible to set different temperatures in each zone. The zone barriers reduce the amount of heat transfer between zones, increasing overall control and enabling the operator to create a temperature gradient across the heated length e.g., zone 1: 500°C, zone 2: 450°C, zone 3: 400°C.

### **WHICH APPLICATIONS REQUIRE A GRADIENT TUBE FURNACE?**

There are many applications that require a gradient tube furnace, such as chemical vapour deposition (CVD), in which a substance is vaporised in the hottest zone, and transported along a work tube via a carrier gas, then condensed upon a substrate in the coolest zone. Gradient tube furnaces are ideal for this process, as each end of the furnace can operate at a different temperature, with a gradient of temperature between them.

[www.carbolite.com/tg](http://www.carbolite.com/tg)