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- | Analysis software which can be used in fully automatic or manual modes for coal ash samples and manual only for biomass and SRF samples.
- | Software zoom function to enable accurate post-test analysis of individual samples with improved resolution
- | One configurable grid assigned to each test piece
- | Temperature controller program set up within the software
- | Space saving embedded computer with Windows IoT Enterprise software runs future proof firmware
- | Default software settings and individual analysis form for coal ash, biomass and SRF
- | An optional work tube integrated lighting system when testing low initial deformation temperature of biomass or SRF samples
- | Lightweight insulation allows quick cooling permitting multiple tests to be completed during the day
- | Automated digital image capture of samples. The frequency of images recorded is set by customer preference, from every 1 °C increment to every 20 °C. The maximum interval for auto analysis is 5°C.


TEST FURNACE FOR ASH FUSIBILITY - CAF G5

TECHNISCHE DETAILS

View inside

1. 1600 °C tube furnace with integral SiC elements
2. External link to embedded PC & software
3. Flow meters for oxidising, reducing gas flow
(dependant on the requirements of the standards)

4. 79 mm inner diameter work tube allows more than 6 samples
5. Digital camera for fast and accurate image recording
6. Gas tight seal for efficient use of gases & safety of operator
7. Automatic temperature programmer with multiple PID control
8. Gas inlets for reducing, oxidising & purge gasses
9. Oxidising or reducing gas selection switch
10. Work tube integrated light for use when testing low 'initial deformation' temperature of biomass and SRF samples (optional)

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TECHNISCHE DETAILS (MODELLEN)

CAF G5

Temperatuursbereik	Up to 1600 °C (1600 °C required for some biomass samples)
Temperature Precision	± 3 °C above 800 °C
Temperature Ramp Rate	7 °C per minute
Temperatuurscontrole	Digital multiple PID terms with gain scheduling and multi offset parameters
Temperature Display	°C
Work Tube dimensions	79 mm internal diameter
Tube material	Mullite
Heating Elements	Silicon carbide x 6
Maximum Sample Load, Manual Analysis	8
Maximum Sample Load, Automatic Analysis	6
Conforms to Standards	BS ISO 540:2008; ASTM D 1857 / D1857M -18); DIN 51730:2007-09; DD CEN/TS 15370-1:2006; PD CEN/TR 15404:2010
Ash Fusibility Determination	Automatic or Manual (Coal & coke: DT, ST, H, FT) Manual only (Biomass / SRF: IST, DT, HT, FT)
Analysetijd	3 runs per working day (including cool down times)
Image Collection	Digital - up to 1 frame per 1 °C rise in temperature
Image Resolution	1280 x 1024 pixels
Gas Requirements: Purge	N2 or CO2
Gas Requirements: Oxidising	CO2 or Air
Gas Requirements: Reducing	CO + CO2 or H2 + CO2
Ventilation	Forced air ventilation
Exhaust	Pipe to be vented into a separate fume hood
Safety	Fail safe gas system and CO alarm supplied
Physical Dimensions (mm)	790 (h) x 505 (w) x 765 (case depth) x 970 (overall depth)
Weight (kg) (furnace)	84
Stroomvoorziening	380 - 415 V, 50/60 Hz two phase 25 A/phase or 220 - 240 V, 50/60 Hz single phase 50 A
Power switching	Solid state relays
Maximum power consumption (kW)	7

Environment Conditions - Operating Conditions	5 °C - 40 °C
Environment Conditions - Relative Humidity	maximum 80 % up to 31 °C decreasing linearly to 50 % at 40 °C
Overtemperature protection	Digital with single high alarm relay

www.carbolite.com/cafg5