TG GRADIENT BUISOVEN REEKS

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HIGH-QUALITY THERMAL INSULATION


LOW ENERGY CONSUMPTION

Carbolite's new tube furnaces are very energy-efficient. Their high-quality thermal insulation is designed for longevity, low case temperature and most of all outstandingly low energy consumption.

Using a Carbolite tube furnace reduces your energy cost as well as the carbon footprint.

UP TO 50% LESS ENERGY REQUIRED!

COMPARISON WITH EQUIVALENT COMPETITOR PRODUCT

Energy consumption for exemplary Carbolite TS1 12/60/600 model (red) and equivalent competitor product (yellow); measured at each temperature after 2 hours soak.

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
A sample could be moved between two temperature zones to achieve a desired heat cycle without waiting for a single zone furnace to heat or cool.

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



TS SCHARNIERENDE BUISOVEN REEKS

THAT'S COOL! LOW CASE TEMPERATURE

Carbolite tube furnaces are designed with user safety in mind.

The robust construction and high-quality thermal insulation ensure that external case temperatures are much lower than other models.

Not only does this help mitigate the risk of injury to the operator, the reduction in the amount of heat escaping the furnace ensures a comfortable working environment, and also means that less energy is wasted during use. When the heat is on, Carbolite can help you keep a cool head!

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TG GRADIENT TUBE FURNACE RANGE - FAQ

WHAT IS A GRADIENT TUBE FURNACE?

A gradient tube furnace is a split tube furnace with multiple heated zones, fitted with zone barriers. The zone

barriers are pieces of insulation that fit between each heating element module, and help reduce the amount of heat transfer between heated zones.

HOW ARE THE GRADIENTS IN A GRADIENT TUBE FURNACE CONTROLLED?

Each heated zone is controlled by its own independent temperature controller, so it is possible to set different temperatures in each zone. The zone barriers reduce the amount of heat transfer between zones, increasing overall control and enabling the operator to create a temperature gradient across the heated length e.g., zone 1: 500°C, zone 2: 450°C, zone 3: 400°C.

WHICH APPLICATIONS REQUIRE A GRADIENT TUBE FURNACE?

There are many applications that require a gradient tube furnace, such as chemical vapour deposition (CVD), in which a substance is vaporised in the hottest zone, and transported along a work tube via a carrier gas, then condensed upon a substrate in the coolest zone. Gradient tube furnaces are ideal for this process, as each end of the furnace can operate at a different temperature, with a gradient of temperature between them.

www.carbolite.com/tg